Aesculap Power Systems

GB

USA

D

F

E

I

P

NL

S

RUS

Care advisory
microspeed® uni motors and motor cables

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GB

USA

microspeed® uni motors and motor cables

Intended use

The microspeed® uni motors are used with a microspeed® uni motor cable connected to microspeed® uni control unit GD670 or microspeed® arthro control unit GD665. The microspeed® uni motors are high-precision, high-performance power units. Consequently, careful processing and care are essential for achieving a long service life.

Overview

| Designation | Art. no. |
|--|----------|
| microspeed® uni motor cable for foot control | GD672 |
| microspeed® uni motor cable with hand control | GD673 |
| microspeed® uni mini 100 motor, mini-Line coupling | GD674 |
| microspeed® uni Hi 150 motor, Hi-Line coupling | GD676 |
| microspeed® uni Hi 100 motor,Hi-Line coupling | GD677 |
| microspeed® uni micro 150 motor, micro-Line coupling | GD678 |
| microspeed® uni micro 100 motor, micro-Line coupling | GD679 |

Safe handling and preparation

- CAUTION**
- Federal law restricts this device to sale by or on order of a physician!**
- Read, follow and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Clean the new product either manually or mechanically prior to the initial sterilization.
 - Store any new or unused products in a dry, clean and safe place.
 - Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.
 - Do not use products that are in open or damaged sterile packaging.

Safe operation

Note
This care advisory does not replace the instructions for use of the microspeed® uni motor system. For operating the microspeed® uni motors and motor cables, always observe the instructions for use of the microspeed® uni motor system (TA022089) or of the microspeed® arthro motor system (TA022156).

Processing

Note
Observe all relevant national regulations and standards with regard to processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note
Up-to-date information on reprocessing can be found on the Aesculap Extra-net at www.aesculap-extra.net.

- Preparations**
- Prior to first mechanical cleaning/disinfection: Mount ECCOS® storage device GB578R (for GD674, GD676, GD677, GD678, GD679) and GB689R (for GD672, GD673) in a suitable basket tray (e.g. JF214R).
 - Mount the storage devices in such a way that the cables can be coiled in a circle with a bending radius of not less than 10 cm. Do not kink the cables.
 - Separate the motor from the motor cable prior to reprocessing.
 - Carry out non-fixating/NaCl-free pre-cleaning immediately after use.
 - Remove all attached components such as adapters, tubes, etc.
 - Reprocess the product immediately after use.
 - Do not clean the product by ultrasound treatment under any circumstances.
 - Put away the product in dry condition.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Only use cleaning and disinfecting agents that do not affect softeners (e.g. silicone), according to the manufacturers' instructions for use.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions regarding concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Preferably use thermal disinfecting processes.

- Mechanical cleaning/disinfecting**
- Pre-cleaning
- Thoroughly pre-rinse under running water.

- Cleaning
- Place the product in its proper position in the ECCOS® storage device immediately after use.
 - Process the product in a single-chamber cleaning/disinfecting machine without ultrasound treatment. Follow the instructions provided by the manufacturer of the machine.
 - Carry out the processing cycle:
 - Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 55 °C.
 - Wash the product for at least 3 min.
 - Neutralize if necessary.
 - Carry out intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: Rinse with distilled, demineralized or fully desalinated water for 5 min, at 93 °C.
 - Complete the program with a drying phase of at least 10 min at a temperature of not more than 120 °C.
 - After completion of the mechanical cleaning/disinfecting cycle, inspect for visible debris all surfaces and places that are difficult to access.
 - Clean manually if necessary.

- Manual cleaning/disinfecting**
- Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
 - Do not clean the product by ultrasound treatment under any circumstances.
 - Do not immerse the product in any fluid. Allow any fluid that has penetrated the device to run out immediately.
 - Clean all surfaces with a soft nylon brush and a suitable cleaning/disinfecting agent.
 - Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
 - Thoroughly brush and rinse surfaces under running water.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.

- Inspect for visible debris all surfaces and places that are difficult to access. If necessary, repeat the cleaning/disinfection process.
- Dry the product with a lint-free cloth.
- Make certain that lumens and channels are dried, too.

- Control, care and inspection**
- Allow the product to cool down to room temperature.
 - Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged and does not show atypical running noises, overheating or excessive vibrations.
 - Set aside the product if it is damaged.
 - After cleaning and before sterilization, always lubricate the motor with oil spray STERILIT® HI GA536:
 - Fully insert the spray can adapter on the nipple at the rear end of the motor.
 - Hold the spray can in a vertical position and push down the spray nozzle valve until oil extrudes at the front end of the handpiece coupling of the motor.
 - Wipe away any excess oil.

- Packaging**
- Always follow the instructions for use of the respective storage devices (e.g. instructions for use TA009721 for the Aesculap ECCOS® storage system).
 - Insert the product in its proper position in the ECCOS® storage device, or put it on a tray in such a way that the product is protected against damage. Observe the weight limit for each tray/container.
 - Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
 - Make certain that the product is wrapped in such a way that recontamination of the product in the period between reprocessing and reuse will be prevented.

Sterilization method and parameters

CAUTION

Damage to motor cables GD672 and GD673 caused by sterilization in sterilizers of small volume! Hot wall surfaces (>150 °C) in the sterilizer can damage the motor cable!

- Protect the motor cable (e.g. by wrapping it in a cloth).

- Sterilize with steam, taking note of the following:
The sterilization has to be carried out according to a validated steam sterilization procedure (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993 and validated according to EN 554/ISO 13683). In the fractioned vacuum process, sterilization must be carried out using the 134 °C/2 bar program with a minimum holding time of 5 minutes.

Note
Steam sterilization at 134 °C for a holding time of 18 min can adversely affect the service life of the product.

- Make certain that the motor cables are not kinked or coiled up too tightly. The coil radius has to be larger than 10 cm.
- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

- Sterilization for the US market**
- Sterilization of the device may be accomplished by steam.
 - Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
 - Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.
- The recommended sterilization parameters are as follows:

| Sterilization method | Temp. | Minimum exposure time | |
|----------------------|------------|-----------------------|-------------------------------|
| | | Wrapped | In a sterile container system |
| Pre-vacuum | 270—275 °F | 4 min | 4 min |

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

- Storage**
- Store processed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

Maintenance

To ensure reliable operation, Aesculap recommends that maintenance be performed after 300 reprocessing cycles, at least once per year.
For services to this end, please contact your national B. Braun/Aesculap agency, see Technical Service.

Technical Service

For service, maintenance or repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints
Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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